

# BMN I-Type

## Set for iliac crest bone marrow aspiration

The BMN I-Type needle is an aspiration needle that is used for bone marrow aspiration of the iliac crest.



- The tip of the needle and the stylet are designed **to allow easier penetration of the bone cortex**
- **Removable control** for variable length
- **The ergonomic handle** enables a better grip and improved manoeuvrability

Part number	Diameter	Length
I 1407M	14 G	9 - 20 - 53 - 70 mm
I 1408M	14 G	19 - 30 - 63 - 80 mm
I 1507M	15 G	9 - 20 - 53 - 70 mm
I 1510M	15 G	39 - 50 - 83 - 100 mm
Included in delivery per package: 10 pieces		