

EasyCut

is a disposable device designed for soft tissues biopsy

The needle is featured with an atraumatic-sharpened cannula made with stainless steel and a “Menghini” tip. The stylet has a pyramidal point. The device is equipped with a special system for vacuum activation by pressing the button on the syringe; therefore it can be used with one hand only.



EasyCut

- The **sliding stopper** of the cannula allows to easily settle the depth to be reached
- The needle includes a **rotating luer lock** in order to give a safe attachment to the syringe
- The device is provided with **inner stylet** that eases the device penetration without skin incision

Ordering Information:

Code	Gauge	Diameter mm	Length
PC01620	16 G	1.60	200 mm
PC01710	17 G	1.40	100 mm
PC01720	17 G	1.40	200 mm
Packaging unit: 10 pieces			

For other sizes, please contact our product specialist