

Ecojekt

Aspiration needle with modified Chiba tip

The Ecojekt aspiration needle has a removable stylet that is perfectly coupled with the needle. The needle's depth stop and depth graduation allow for safe control of the penetration depth. The internal echoic marking prevents cancer cells from spreading.



- **Removable mandrin**
- **Internal echoic marker** prevents cancer cells from spreading
- **Safe control of the penetration** by means of depth marker and stopper

Part number	Gauge	Length
E1715	17 G	150 mm
E1815	18 G	150 mm
E1825	18 G	250 mm
E1910	19 G	100 mm
E1915	19 G	150 mm
E1918	19 G	180 mm
E2010	20 G	100 mm
E2015	20 G	150 mm
E2020	20 G	200 mm
E2110	21 G	100 mm
Included in delivery per package: 20 pieces		

Part number	Gauge	Length
E2115	21 G	150 mm
E2120	21 G	200 mm
E2130	21 G	300 mm
E2205	22 G	50 mm
E22075	22 G	75 mm
E2210	22 G	100 mm
E2215	22 G	150 mm
E2220	22 G	200 mm
E2315	23 G	150 mm
Included in delivery per package: 20 pieces		