


SemiCut


Disposable device designed for soft tissues biopsy such as breast, lung, liver, and kidney

After cannula loading, SemiCut can be used with one hand, leaving the other free to handle over devices such as ultrasound guide etc. The stylet is equipped with a sample notch of 20 mm and an ultra-sharp triple-face tip to penetrate easily also fibrous tissues. The echogenic marker on the cannula surface allows easy determination of the depth of insertion.



- Biopsy specimen can be taken **with two different lengths: 10 mm or 20 mm**
- **Lightest semiautomatic device on the market**

SemiCut			
Partnumber	Needle Gauge	Diameter mm	Needle Length L
PDO1406	14 G	2.1	60 mm
PDO1410	14 G	2.1	100 mm
PDO1420	14 G	2.1	200 mm
PDO1606	16 G	1.6	60 mm
PDO1610	16 G	1.6	100 mm
PDO1615	16 G	1.6	150 mm
PDO1620	16 G	1.6	200 mm
PDO1806	18 G	1.2	60 mm
PDO1810	18 G	1.2	100 mm
PDO1815	18 G	1.2	150 mm
PDO2007	20 G	0.9	70 mm
PDO2010	20 G	0.9	100 mm
PDO2015	20 G	0.9	150 mm
Included in delivery per package: 10 pieces			



SemiCut & SemiCut MRI with Coaxial				
Partnumber	Needle Gauge	Needle Length L	Coaxiale Gauge	Coaxial Length l
PDOPP1410	14 G	100 mm	13 G	82 mm
PDOPP1415	14 G	150 mm	13 G	132 mm
PDOPP1420	14 G	200 mm	13 G	182 mm
PDOPP1610	16 G	100 mm	15 G	83 mm
PDOPP1615	16 G	150 mm	15 G	133 mm
PDOPP1620	16 G	200 mm	15 G	183 mm
PDOPP1810	18 G	100 mm	17 G	83 mm
PDOPP1815	18 G	150 mm	17 G	133 mm
PDOPP1820	18 G	200 mm	17 G	183 mm
PDOPP2010	20 G	100 mm	19 G	84 mm
PDOPP2015	20 G	150 mm	19 G	134 mm
PDOPP2020	20 G	200 mm	19 G	184 mm
PDOPPMR1410	14G	100 mm	13 G	82 mm
PDOPPMR1415	14G	150 mm	13 G	132 mm
PDOPPMR1610	16G	100 mm	15 G	83 mm
PDOPPMR1615	16G	150 mm	15 G	133 mm
Included in delivery per package: 10 pieces				



For other sizes, please contact our product specialist