


Top Cut /Top Cut M

Tru-Cut Nadel for automated biopsy systems

The TOP CUT and the TOP CUT M is a needle that has a Menghini Tru-Cut tip. The needle has a depth stopper that makes it easier to position the needle within the automatic system.

- **Menghini-Tru-Cut tip** with marker
- **22 mm penetration / 11mm penetration**



Top Cut M				
Part number	Gauge	Diameter mm	Length	System compatibility
TOP 1410 M	14G	2.10	100 mm	Bard® Magnum® biopsy gun
TOP 1413 M	14G	2.10	130 mm	Bard® Magnum® biopsy gun
TOP 1616 M	16G	1.60	160 mm	Bard® Magnum® biopsy gun
TOP 1813 M	18G	1.20	130 mm	Bard® Magnum® biopsy gun
TOP 1820 M	18G	1.20	200 mm	Bard® Magnum® biopsy gun
TOP 1825 M	18G	1.20	250 mm	Bard® Magnum® biopsy gun
Included in delivery per package: 10 pieces				

Top Cut				
Part number	Gauge	Diameter mm	Length	System compatibility
TOP 1610	16G	1.60	100 mm	Pro-Mag™ 2.2 biopsy gun
TOP 1615	16G	1.60	150 mm	Pro-Mag™ 2.2 biopsy gun
TOP 1810	18G	1.20	100 mm	Pro-Mag™ 2.2 biopsy gun
TOP 1820	18G	1.20	200 mm	Pro-Mag™ 2.2 biopsy gun
TOP 1825	18G	1.20	250 mm	Pro-Mag™ 2.2 biopsy gun
Included in delivery per package: 10 pieces				

For other sizes, please contact our product specialist