

TRAPSYSTEM

For bone marrow biopsy

The TRAPSYSTEM device makes it possible to perform bone marrow biopsy without using the conventional deflection method for obtaining a bone marrow sample. Thanks to its specially developed cannula, it catches the sample in its full size and undamaged when inserted into the "J" needle.

- **"Fish mouth"-like tip**
- **"J" needle with** outer centimetre graduation
- **Stylet with pyramid-shaped tip**
- **Extraction cannula with centimetre graduation** and dull tip
- **Ergonomic handle**
- **Safe** sampling
- **Easy maneuverability**
- **Fast** method



Part number	Gauge	Length
TRAPJ 0810	8 G	100 mm
TRAPJ 0815	8 G	150 mm
TRAPJ 1110	11 G	100 mm
TRAPJ 1115	11 G	150 mm
TRAPJ 1306	13 G	60 mm
TRAPJ 1310	13 G	100 mm
Included in delivery per package: 10 pieces		