

Everlift®

Submucosal lifting agent

EverLift Submucosal Lifting Agent is indicated for use in gastrointestinal endoscopic procedures for submucosal lift of polyps, adenomas, early-stage cancer, or other gastrointestinal lesions prior to excision with a snare or other appropriate endoscopic device.

EverLift is cellulose-based, which means the solution stays in place longer as it is not absorbed as quickly as other agents.

EverLift stains the submucosa with a deep blue colour, providing a strong contrast to the white colour of the muscle layer.



- Pre-filled syringe means **no preparation time**
- Cellulose based & hypertonic, allowing **2 mm lift for 60 minutes**
- **The only pre-filled submucosal lifting agent on the market, available in both 5 ml and 10 ml sizes**
- The syringe can be connected directly to all endoscopic injection needles via the Luer lock connection. **23G needle recommended.**

Article number	GIS-55	GIS-59
Content per syringe	5 ml	10 ml
Scope of delivery per package: 10 units		

