


3-level balloon catheter

For stone extraction

The 3-level balloon catheter is used for stone extraction. It consists of one latex-free balloon. The catheter is triple-lumen and 200 cm long.



- **Latex-free** balloon
- **Improved X-ray contrast**
- **Triple-lumen catheter** 200 cm long
- **Can be used with guidewire** .035 inch
- **Including three corresponding** Filling syringes

Part number	Design	Catheter diameter	Catheter length	Balloon diameter when filled		
N135030PRO	Triple-lumen	7 Fr	200 cm	9 mm	13 mm	16 mm
Included in delivery per package: 1 piece						