


# Biopsy forceps

## Coated, for single use

Coated biopsy forceps for single use. For working channels 2.8 mm in diameter and larger. Available with or without pins.



- **Coated**
- For working channels **2.8 mm in diameter and larger**
- 2.3 mm  $\varnothing$  incl. handle
- **Available in two versions:**
  - Without pins
  - With pins

	Without pins		With pins	
<b>Part number</b>	<b>NEX0222-B</b>	<b>NEX0222-G</b>	<b>NEX0222D-B</b>	<b>NEX0222D-G</b>
<b>Diameter <math>\varnothing</math></b>	<b>2.3 mm</b>	<b>2.3 mm</b>	<b>2.3 mm</b>	<b>2.3 mm</b>
<b>Length</b>	<b>160 mm</b>	<b>230 mm</b>	<b>160 mm</b>	<b>230 mm</b>
<b>Included in delivery per package: 10 pieces</b>				