

Cysto-Gastro Kit

For the puncture of pancreatic pseudocysts

By using endoscopic ultrasound or fluoroscopy pancreatic pseudocysts can be punctured transgastrally or transduodenally with the cysto-gastro set. The set consists of a needle and outer tube with a diathermic tip, two HF connectors for the needle papillotome.

The Cystostome Fr. 6 is specially designed for use in smaller working channels. This consists of an outer tube with a corresponding diathermic metal tip and an attachment for a .035 inch guidewire.

Following the puncture of the pseudocysts, drainage stents can be installed for all of the sets in the stomach or the duodenum using the guidewire.



- consisting of **needle and outer tube**
- **two HF connectors**
- **drain** installed using guidewire

Part number	FR	Length
E30006 345	6	180 cm
E30008 341	8.5	180 cm
E30010 340	10	180 cm
Included in delivery per package: 1 piece		