

EasyCollect

Oval collection bag for single use

EasyCollect for single use. The collection bag can be turned in all directions and is latex-free.

- **Ovaler collection bag**
- **Can be turned** in all directions
- **Latex-free**
- **Ergonomic** plastic handle
- **Individual, sterile packaging**
- **Available in two sizes**
 - Collection bag 15 mm
 - Collection bag 35 mm



Part number	NE1518-B15	NE1522-G35
Bag size	15 mm	35 mm
Outer diameter	1.8 mm	2.3 mm
Length	160 cm	230 cm
Included in delivery per package: 10 pieces		