

Hemoclip

For haemostasis

The Hemoclip is suitable for application with a flexible endoscope. Primarily for stopping bleeding of various causes in the upper and lower digestive tract and for postoperative wound suture.

- **Repeated** opening and closing possible
- **360° bi-directional synchronous** rotational control
- **Excellent** clamping strength
- **Unique opening span with optional shifting function between:**
 - 11 mm
 - 16 mm
- **Removable function:** Clip can be pro-actively removed from the tissue in special cases



Part number	Length	Clip angle	Breadth	Diameter
AG-51044-2300-090-16	230 cm	90°	16 mm	2.6 mm
AG-51044-2300-135-16	230 cm	135°	16 mm	2.6 mm
AG-5104-2300-90	230 cm	90°	11 mm	2.6 mm
AG-5104-2300-135	230 cm	135°	11 mm	2.6 mm
AG-55504-2300-090-11-2*	230 cm	90°	11 mm	2.6 mm
AG-55504-2300-135-11-2*	230 cm	135°	11 mm	2.6 mm
AG-55504-2300-090-16-2*	230 cm	90°	16 mm	2.6 mm
AG-55504-2300-135-16-2*	230 cm	135°	16 mm	2.6 mm
Included in delivery per package: 5 pieces				

* Rechargeable, including 2 extra clips