


# Lithotripter Set

## Reusable

Mechanical lithotripter set for a combination of planned and emergency lithotripsy in one device. The attachment of the spiral and the securing of the ratchet ensure ease of use and the set is compatible with all commonly used Endo-Flex lithotripter baskets. The set is available in two different options and reusable.



- **Easy to use**
- Emergency and planned lithotripsy **combined in one device**
- **Compatible** with lithotripter baskets

Part number	Description	Use
120610	SET 1: mechanical lithotripter set type 3	Reusable
120611	SET 2: lithotripter handle incl.2x lithotripter spirals	Reusable
120600	Lithotripter handle type 3	Reusable
E120602	Lithotripter spiral Ø 2.6 mm, L 195 cm	Single use
E120604	Lithotripter spiral Ø 4.0 mm, L 100 cm, for type 3	Single use
E122***N	Lithotripter basket, diamond shaped, 4 multifile wires, Ø 2.6 mm, L 400 cm, basket 50 mm	Single use
E124***N	Lithotripter basket, twisted, 6 multifile wires, Ø 2.6 mm, L 400 cm, basket 50 mm	Single use
E128336N	Lithotripter basket, diamond shaped, 6 multifile wires, Ø 3.4 mm, L 400 cm, basket 60 mm	Single use

Included in delivery per package: 1 piece

\*\*\* For other sizes, please contact our product specialist