

Multiband Ligature System

“MBL”

The Multiband Ligature System can be used with all endoscopes up to 180 cm in length. Very easy introduction of the lasso into the working channel saves you a lot of preparation time.



- Preloaded for **6 ligatures**
- **Very easy introduction of the lasso** into the working channel
- **Side LL rinse attachment** with LL cap
- For working channels **2.8 mm in diameter and larger**
- Can be used with **all endoscopes up to 180 cm in length**
- **Individually** packaged

Part number	MBL-01-06
Included in delivery per package:	1 piece