


Tuloc® / Tuloc® Premium

Correctable localisation system

The proven and reliable Tuloc® localisation system is used for the preoperative marking of non-palpable, suspicious breast lesions.

The system is repositionable and makes the intraoperative localisation of findings easier by means of a depth marking on the wire. Tuloc® Premium also features an ergonomic handle for safe, one-handed operation.

For use in MRI we recommend the MRI Tuloc® exclusively.



- **Repositionable marking wire**
- **Monofilament wire** for best possible form and pressure stability
- **Special material and loop shape** ensure good visibility in ultrasound and palpability during surgery
- **Fixation piece** allows for safe preoperative fixation of the wire

Part number	Description	Gauge	Diameter	Length
271661	Tuloc	20	0.95 mm	90 mm
271663	Tuloc	20	0.95 mm	120 mm
271881	Tuloc Premium	20	0.95 mm	90 mm
271883	Tuloc Premium	20	0.95 mm	120 mm
601649*	MR Tuloc	20	0.95 mm	90 mm
601651*	MR Tuloc	20	0.95 mm	120 mm

Units per package: 10 pieces

*MRI Tuloc® is not stocked at Zuchwil, but can be special ordered from the supplier