

KliniTray™ vulva plate

Transportation and X-ray system for the surgical excision of vulval cancer

The KliniTray™ vulva plate tissue transportation and X-ray system is an innovative tissue fixation system used to transport, examine, document and prepare vulvar samples following the surgical excision of vulval cancer.



With the help of the KliniTray™ vulva plate, the tissue can be transported safely in a determined three-dimensional direction and examined using radiological equipment and ultrasound.

The system also comes with a light plastic tray printed with a visual depiction of the vulva and the directional markings “dorsal” and “ventral”, which map out the location of the specimen in the patient’s body from the surgeon’s perspective.

- **Safe anatomical allocation** of the resection margins
- **Safe, air- and water-tight** transport container
- **Lightweight plastic tray** for tissue fixation
- Clear localization of findings thanks to **radiolucent directional markers**

Part number	Description	Dimensions
2001300	KliniTray™ vulva-plate, sterile, with 5 metal pins	18 x 20 x 1 cm
1000350	Plastic container for KliniTray™ vulva plate, non-sterile	3000 ml
Included in delivery per package: 20 pieces		