

Chiba needle

With high level of echogenicity

The Chiba needle is produced with a special cut that guarantees a high level of echogenicity.



- **High level of echogenicity** from a special cut
- Mostly **atraumatic puncture**

Part number	Diameter gauge / mm	Length	Quantity per package
200S 095150	19.5 G / 0.95 mm	150 mm	10 pieces
200S 095220	19.5 G / 0.95 mm	220 mm	10 pieces
200S 070090	22 G / 0.70 mm	90 mm	25 pieces
200S 070120	22 G / 0.70 mm	120 mm	10 pieces
200S 070150	22 G / 0.70 mm	150 mm	10 pieces
200S 070220	22 G / 0.70 mm	220 mm	10 pieces